

Forward Transfection Using Lipofectamine® RNAiMAX

Use this procedure to forward transfect Stealth™ RNAi or siRNA into mammalian cells in a 6- well format (for other formats, see Scaling Up or Down Transfections Table) . In forward transfections, cells are plated in the wells, and the transfection mix is generally prepared and added the next day. All amounts and volumes are given on a per well basis.

Reagents:

Lipofectamine® RNAiMAX (Thermo Fisher Scientific # 13778150); RNAi duplex solution; Sterile tubes;

Procedures:

1. One day before transfection, plate cells in 2 ml of growth medium without antibiotics such that they will be 50-60% confluent at the time of transfection.
2. For each well to be transfected, prepare RNAi duplex-Lipofectamine® RNAiMAX complexes as follows:
 - Dilute 20 pmol RNAi duplex in 50 µl Opti-MEM® I Reduced Serum Medium without serum. Mix gently.
 - Mix Lipofectamine™ RNAiMAX gently before use, then dilute 1 µl in 50 µl Opti-MEM® I Reduced Serum Medium. Mix gently.
 - Combine the diluted RNAi duplex with the diluted Lipofectamine™ RNAiMAX. Mix gently and incubate for 10 minutes (no more than) at room temperature.
3. Add the RNAi duplex-Lipofectamine® RNAiMAX complexes (100 µl) to each well containing cells. This gives a final RNA concentration of ~10 nM per well. Mix gently by rocking the plate back and forth.
4. Incubate the cells 48-72 hours at 37°C in a CO₂ incubator until you are ready to assay for gene knockdown. Medium may be changed after 4-6 hours.

Table: Scaling Up or Down Transfections for Various Formates

<i>Components</i>	10-cm	6-cm	6-well	12-well	24-well	96-well	Note
<i>Lipofectamine® RNAiMAX (μl)</i>	5 μl	2 μl	1 μl	0.4 μl	0.2 μl	0.05 μl	
<i>Opti-MEM® I Medium for Lipo dilution</i>	200 μl	100 μl	50 μl	20 μl	10 μl	2.5 μl	**1
<i>RNAi duplex (20 μM stock)</i>	5 μl	2 μl	1 μl	0.4 μl	0.2 μl	0.05 μl	**2
<i>Opti-MEM® I Medium for siRNA dilution</i>	200 μl	100 μl	50 μl	20 μl	10 μl	2.5 μl	
<i>Final siRNA concentration (nM)</i>	10 nM	10 nM	10 nM	10 nM	10 nM	10 nM	

**1: Other medium without serum/antibiotic can be used to replace Opti-MEM.

**2: The final concentration of siRNA is 10 nM in this table. The final concentration can be reduced to 1 nM without reduction in knockdown efficacy for most of the siRNA we have tested. A titration test (for example, 1, 5, 10 nM) is recommended for new siRNA. The minimal amount with good knockdown efficiency is recommended.

Modified from: <https://www.thermofisher.com/no/en/home/references/protocols/cell-culture/transfection-protocol/rnaimax-forward-transfections-lipofectamine.html> and https://tools.thermofisher.com/content/sfs/manuals/Lipofectamine_RNAiMAX_Reag_protocol.pdf

Edited by Yang Jin 17/11/2017