

## Iron oxidoreductase assay

### 1. Reaction Buffer

	Final conc. (mM)
MES	25
MOPS	25
NaCl	140
Glucose	5
KCl	5.4
CaCl <sub>2</sub>	1.8
MgCl <sub>2</sub>	0.8
Adjust pH to 7.4	

###MES and MOPS can be replaced with 16,7 mM HEPES.

### 2. Fe(III)-NTA3 (*iron(III)-nitrilotriacetate (NTA)*) stock solution (1 mM Fe(III) plus 3 mM NTA)

\*\*\*Ferric nitrate (FeNO<sub>3</sub>) or FeCl<sub>3</sub> was dissolved in 100 mM HCl. Disodium nitrilotriacetate (Na<sub>2</sub>NTA) was dissolved in H<sub>2</sub>O (If Nitrilotriacetic acid is used, 10 M NaOH solution is needed to help NTA acid dissolution). The respective solutions were mixed to achieve a molar ratio of 1:3 of Fe-NTA. The pH was adjusted to 7.4 with sodium bicarbonate with constant stirring (I prepared Fe-NTA stock solution at 10 mM). (Ref. *Awai M, Narasaki M, Yamanoi Y, Seno S. Induction of diabetes in animals by parenteral administration of ferric nitrilotriacetate. A model of experimental hemochromatosis. Am J Pathol 1979;95:663–72.*)

### 3. Ferrozine

Desired amount of ferrozine can be dissolved in either H<sub>2</sub>O or Reaction Buffer (without Fe-NTA<sub>3</sub>) just before use. Keep the solution in dark.

#### Procedures:

1. Transfect cells (in 6-well or 12-well plate) with target DNA
2. Forty-eight hours after transfection, cells are washed once with phosphate-buffered saline (PBS), pH 7.2.
3. Prepare Iron Uptake Buffer (IUB, containing 50-100  $\mu\text{M}$  Fe-(NTA)<sub>3</sub> and 1 mM (at least 200  $\mu\text{M}$ ) ferrozine) by combine Reaction buffer, ferrozine and Fe-NTA<sub>3</sub>.
4. Incubate cells IN DARK in Iron Uptake Buffer (800  $\mu\text{l}$ /well of 12-well plate) 37 °C for **20-60 min**. (\*\*\*) Once ferrozine is mixed with Fe-NTA and buffer, the mixture must be kept from light.)
5. After incubation, the Fe<sup>2+</sup>-ferrozine complex was detected by measuring the absorbance at 562 nm and converted to Fe<sup>2+</sup> generated by using the extinction coefficient 27.9 mM<sup>-1</sup>cm<sup>-1</sup>.

Cell number : Normally an 80% confluent well of HEK293T cells in a 6 or 12 well plate should give you a reliable assay after 1 hr. incubation.

\*\*\*\* Original protocol suggests that all reagents should be prepared before use. In my experience, **Fe(III)-NTA<sub>3</sub> can be stored in fridge for several months. Reaction buffer (without Fe(III)-NTA<sub>3</sub> and ferrozine can be kept at -20 degree for very long time. Normally 4x reaction buffer is prepared to minimize solution volume for storage. )**